

## Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella O 27</b>
Catalogue No. (Package size)	TR 1313 (1 ml)
Lot No.	<b>14 01 14</b>
Expiry date	<b>2018-08</b>
Intended use	Enteroclon Anti-Salmonella O 27 records all Salmonella serovars of group B and D <sub>3</sub> as long as they contain the O 27 antigen.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

22.01.14

  
\_\_\_\_\_  
Head of Quality Control