

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 28
Catalogue No. (Package size)	TR 1336 (1 ml)
Lot No.	01 01 12
Expiry date	2017-01
Intended use	Enteroclon Anti-Salmonella O 28 records all Salmonella strains, as long as they contain antigen O 28.
Specificity	No agglutination with heterologous Salmonella O antigens and Salmonella H antigens
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 08.02.12


Head of Quality Control