

### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Enteroclon Anti-Salmonella O 38</b>
Catalogue No. (Package size)	TR 1344 (1 ml)
Lot No.	<b>119 08 13</b>
Expiry date	<b>2018-08</b>
Intended use	Enteroclon Anti-Salmonella O 38 records all Salmonella strains, as long as they contain antigen O 38.
Specificity	No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

16.01.14

  
Head of Quality Control