


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Enteroclon Anti-Salmonella O 42
Catalogue No. (Package size)	TR 1348 (1 ml)
Lot No.	129 08 12
Expiry date	2015-08
Intended use	Enteroclon Anti-Salmonella O 42 records all Salmonella strains, as long as they contain antigen O 42.
Specificity	No other agglutinations with heterologous Salmonella O antigens and Salmonella H antigens are known.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

16.08.12


Head of Quality Control