

### Certificate of Quality Control

|                                 |   |
|---------------------------------|---|
| Manufacturer                    | <b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin,<br>Berliner Allee 317-321, 13088 Berlin, Germany |
| Product                         | <b>Enteroclon Anti-Salmonella O 45</b>  |
| Catalogue No.<br>(Package size) | TR 1351 (1 ml)  |
| Lot No.                         | <b>114 08 12</b>  |
| Expiry date                     | <b>2017-08</b>  |
| Intended use                    | Enteroclon Anti-Salmonella O 45 records all Salmonella strains, as long as they contain antigen O 45.                 |
| Specificity                     | No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.                                   |
| Preservatives                   | Sodium azide 0.9 mg/ ml   |
| Storage                         | 2...8 °C  |

Berlin, 23.08.12

  
\_\_\_\_\_  
Head of Quality Control