

### Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,  
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella O 67**

Catalogue No.  
(Package size) TR 1371 (1 ml)

Lot No. **50 06 10**

Expiry date **2015-06**

Intended use Enteroclon Anti-Salmonella O 67 records all Salmonella strains, as long as they contain antigen O 67.

Specificity No agglutination with heterologous Salmonella O antigens and Salmonella H antigens.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin, 29.07.10

  
Head of Quality Control