

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:c
Catalogue no.	TR 1403
Package size	1 ml
Lot no.	1580219
Expiry date	2022-02
Intended use	Anti-Salmonella H:c records all <i>Salmonella</i> strains, as long as they contain antigen H:c in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,
place, date

27.02.2019

signature Head of Quality Control

