

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:d
Catalogue no.	TR1404
Package size	1 ml
Lot no.	9221121
Expiry date	2026-10
Intended use	Anti-Salmonella H:d records all <i>Salmonella</i> strains, as long as they contain antigen H:d in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 11.01.22
place, date

F. Luch
signature Head of Quality Control