

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Salmonella H:g</b>
Catalogue no.	TR 1406
Package size	1 ml
<b>Lot no.</b>	<b>9911119</b>
<b>Expiry date</b>	<b>2022-11</b>
Intended use	Anti-Salmonella H:g records all <i>Salmonella</i> strains, as long as they contain antigen H:g (e.g. H:f,g; H:f,g,s; H:f,g,t; H:g,m; H:g,m,s; H:g,p; H:g,p,u; H:g,q; H:g,s,t; H:g,t; H:g,z <sub>51</sub> ) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 06.12.19  
place, date

J. Lendner  
signature Head of Quality Control