

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Salmonella H i**

Catalogue No.
(Package size) TR 1410 (1 ml)

Lot No. **103 07 13**

Expiry date **2016-07**

Intended use Enteroclon Anti-Salmonella H i records all Salmonella strains, as
long as they contain antigen H i in slide agglutination.

Specificity No agglutination with heterologous Salmonella H antigens and
Salmonella O antigens.

Preservatives Sodium azide 0.9 mg/ ml

Storage 2...8 °C

Berlin,

05.09.13


Head of Quality Control