

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

Product	Anti-Salmonella H:L
Catalogue no.	TR 1412
Package size	1 ml
Lot no.	8180917
Expiry date	2020-09
Intended use	Anti-Salmonella H:L records all <i>Salmonella</i> strains, as long as they contain antigen H:l (e.g. H:l,v; H:l,w; H:l,z ₁₃ ; H:l,z ₂₈ ; H:l,z ₁₃ ,z ₂₈) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin,

09.11.17

place, date

F. Ludwig

signature Head of Quality Control