

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany

T: +49 30 927 030-0
F: +49 30 927 030-30
info@sifin.de
www.sifin.de

| | |
|--------------------|---|
| Product | Anti-Salmonella H:L |
| Catalogue no. | TR 1412 |
| Package size | 1 ml |
| Lot no. | 1020119 |
| Expiry date | 2022-01 |
| Intended use | Anti-Salmonella H:L records all <i>Salmonella</i> strains, as long as they contain antigen H:ℓ (e.g. H:ℓ,v; H:ℓ,w; H:ℓ,z ₁₃ ; H:ℓ,z ₂₈ ; H:ℓ,z ₁₃ ,z ₂₈) in slide agglutination. |
| Specificity | No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens. |
| Preservatives | Sodium azide 0.9 mg/ml |
| Storage | 2 °C...8 °C |

Berlin,
place, date

27.02.2019

signature Head of Quality Control

