

Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

Product	Anti-Salmonella H:L
Catalogue no.	TR1412
Package size	1 ml
Lot no.	8201020
Expiry date	2023-09
Intended use	Anti-Salmonella H:L records all <i>Salmonella</i> strains, as long as they contain antigen H:ℓ (e.g. H:ℓ,v; H:ℓ,w; H:ℓ,z ₁₃ ; H:ℓ,z ₂₈ ; H:ℓ,z ₁₃ ,z ₂₈) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 04. 11. 2020
place, date

F. Muderig
signature Head of Quality Control