

Certificate of Quality Control

Manufacturer **SIFIN** Institut für Immunpräparate und Nährmedien GmbH Berlin,
Berliner Allee 317-321, 13088 Berlin, Germany

Product **Enteroclon Anti-Shigella sonnei S- and F-form (Phase I
and II)**

Catalogue No. TR 2014 (1 ml)
(Package size)

Lot No. **16 01 14**

Expiry date **2016-02**

Intended use Enteroclon Anti-Shigella sonnei S- and F-form (Phase I and II)
records the strains of Shigella sonnei S- (Phase I) and F-form
(Phase II) in slide agglutination.

Specificity No agglutination with antigens of Shigella flexneri, Shigella
dysenteriae (type 1 to 10) and Shigella boydii (type 1 to 15).

Preservatives Sodium azide 0.9 mg / ml

Storage 2...8 °C

Berlin,

11.02.14


Head of Quality Control