

## Certificate of Quality Control

CERTIFIED QUALITY MANAGEMENT SYSTEM IN  
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

<b>Product</b>	<b>Anti-Salmonella H:m</b>
Catalogue no.	TS1413
Package size	1 ml
<b>Lot no.</b>	<b>7230920</b>
<b>Expiry date</b>	<b>2025-08</b>
Intended use	Anti-Salmonella H:m records all <i>Salmonella</i> strains, as long as they contain antigen H:m (e. g. H:g,m; H:g,m,s,t; H:m,t) in slide agglutination.
Specificity	No agglutination with heterologous <i>Salmonella</i> H antigens and <i>Salmonella</i> O antigens.
Preservatives	Sodium azide 0.9 mg/ml
Storage	2 °C...8 °C

Berlin, 30.09.2020

place, date



signature Head of Quality Control