


Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin
Product	Test Serum Anti-Salmonella H x
Catalogue No. (Package size)	TS 1422 (1 ml)
Lot No.	117 08 12
Expiry date	2017-08
Intended use	Testserum Anti-Salmonella H x records all Salmonella strains, as long as they contain antigen H x (e.g. e,n,x, e,n,x,z ₁₅) in slide agglutination.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin, 14.07.12


Head of Quality Control