

Certificate of Quality Control

Manufacturer	SIFIN Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	Test Serum Anti-Salmonella H z₁₅
Catalogue No. (Package size)	TS 1428 (1 ml)
Lot No.	149 10 123
Expiry date	2018-10
Intended use	Test Serum Anti-Salmonella H z ₁₅ records all Salmonella strains, as long as they contain in an antigen combination the antigen H z ₁₅ (e.g. e ₇ n ₁ x ₁ z ₁₅ ; e ₇ n ₁ z ₁₅) in slide agglutination.
Specificity	No agglutination with heterologous Salmonella H antigens and Salmonella O antigens.
Preservatives	Sodium azide 0.9 mg/ ml
Storage	2...8 °C

Berlin,

21.11.13


Head of Quality Control