


### Certificate of Quality Control

Manufacturer	<b>SIFIN</b> Institut für Immunpräparate und Nährmedien GmbH Berlin, Berliner Allee 317-321, 13088 Berlin, Germany
Product	<b>Control Serum Anti-Yersinia enterocolitica O 9 for the Widal reaction</b>
Catalogue No. (Package size)	TS 1735 (1 ml)
Lot No.	<b>100 07 13</b>
Expiry date	<b>2018-07</b>
Intended use	The control serum agglutinates on the MTP with the <b>coloured</b> test antigen Yersinia enterocolitica O 9 (SIFIN) with a titre of 1:1600. Deviations of $\pm$ one titre step are permitted.
Specificity	The serum control (antigen addition is substituted by physiological saline solution) must not show any agglutination.
Preservatives	Sodium azide 0.9 mg / ml
Storage	2...8 °C

Berlin,

23.07.13

  
Head of Quality Control